



Canadian Centre for
Computational
Genomics


C3G Analysis Workshop: RNA-Seq

GenPipes Files: Creation & Validation

22-23 January 2019


GenPipes Overview

Experimental Processing:



GenPipes Overview

GenPipes Requirements:



GenPipes requires setting up 3 files:

- Configuration “ini” file: contains parameters used by tools in the pipeline
- Readset file: contains details about your samples
- Design file: contains details about sample comparisons. **NOT needed by all pipelines.**



What are 'configuration/ini' files and why do we need them?



- A “configuration” file, also know as an “ini” file due to its file extention, is a file that stores all the parameters needed by tools in the pipelines.
- A pipeline is a multi-step framework. At every step, tools require several parameters to be set. Having to type all the parameters on the command line as you call the tool would be complicated and messy.
- As a solution, parameters are stored in these ini files.
- To open an ini file and check out its content try:

```
cat $MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.base.ini
```

- ini content:

```
[DEFAULT]
assembly_dir =
$MUGQIC_INSTALL_HOME/genomes/C3G_workshop/%(scientific_name)s.%(
assembly)s
scientific_name=Homo_sapiens
assembly=GRCh38.chr19
assembly_synonyms=hg38

module_R=mugqic/R_Bioconductor/3.5.0_3.7

[picard_sam_to_fastq]
cluster_cpu=-N 1 -n 4
cluster_queue=--mem=16G
```

Configurations Files:



- ini files come with each pipeline in:

`$MUGQIC_PIPELINES_HOME/pipelines/<pipeline_name>/<pipeline_name>*.ini`

```
## for the rnaseq pipeline:  
ls $MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq*.ini
```

- ini files for different species are available in: `$MUGQIC_PIPELINES_HOME/resources/genomes/config/`

```
## to check available species ini files:  
ls $MUGQIC_PIPELINES_HOME/resources/genomes/config/
```

Creating your own configuration/ini file

- If you are interested in editing parameters in the ini file, you can create your own ini file that over-writes the set parameters.
- For example, to edit the number of threads in the star_align step in the rnaseq pipeline from 20 to 5, you can create a file, myini.ini, and in it you can change the parameters. You need to indicate the [section].

```
[star_align]  
threads=5
```

How to use the configuration/ini file



- Configuration files are provided to the pipelines through the flag `-c` or `--config`, as follows:

```
rnaseq.py -c $MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.base.ini
```

- For every pipeline, there are several provided ini files. The `<pipeline>.base.ini` contains all the parameters but is optimized for our internal server.
- To run a pipeline on cedar, you need to add the `<pipeline>.cedar.ini`:

```
rnaseq.py -c $MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.base.ini  
$MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.cedar.ini
```

- To use your own ini, you add it at the end:

```
rnaseq.py -c $MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.base.ini  
$MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.cedar.ini myini.ini
```


- To align to mouse mm9, instead of human, you need to add the mouse specific ini:

```
rnaseq.py -c $MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.base.ini  
$MUGQIC_PIPELINES_HOME/pipelines/rnaseq/rnaseq.cedar.ini  
$MUGQIC_PIPELINES_HOME/resources/genomes/config/Mus_musculus.mm9.ini
```

- The parameters in later ini files will over-write those in earlier ones.

GenPipes Overview

GenPipes Requirements:



GenPipes requires setting up 3 files:

- Configuration “ini” file: contains parameters used by tools in the pipeline
- Readset file: contains details about your samples
- Design file: contains details about sample comparisons. **NOT needed by all pipelines.**



What is the 'Readset' file?



- A Readset file is a file that contains information about your samples.
- It is a **tab-separated** file, that contains the following columns:

```
Sample Readset Library RunType Run Lane Adapter1 Adapter2 QualityOffset BED FASTQ1 FASTQ2 BAM
```

- Some fields are optional and others are mandatory ([check here](#)).
- Sample vs Readset:
 - Samples can be thought of as biological replicates while Readsets as technical replicates. If a library was divided across lanes, then each lane would be a “Readset” of the original “Sample”. In most pipelines, readsets belonging to the same sample are merged under a single sample name.

Creating and Using the Readset file

•How to create the Readset file:

- There are several ways to create a readset file. The most basic way is to create it in Excel manually and copy it to the server.
- Another way is to use the csvToReadset.R script if you have a csv file from Genome Quebec


•How to use a Readset file

- a Readset file is provided to the pipeline using the (-r or --readsets)

```
rnaseq.py -r myreadset.txt
```

GenPipes Overview

GenPipes Requirements:



GenPipes requires setting up 3 files:

- Configuration “ini” file: contains parameters used by tools in the pipeline
- Readset file: contains details about your samples
- Design file: contains details about sample comparisons. **NOT needed by all pipelines.**

Design Files



- The Design file is needed by **SOME** of the pipelines.
- It contains data about sample comparisons. The format of the file depends on the pipeline. Consult the pipeline of interest for [more information](#).
- In general, the design file is **tab-separated**, and has the following structure:

RNA-Seq

Sample	Contrast1	Contrast2
sampleA	1	1
sampleB	2	0
sampleC	0	2

ChIP-Seq

Sample	Contrast1,N	Contrast2,B
sampleA	1	1
sampleB	2	0
sampleC	0	2

- The first column is the **Sample** column. Every column after that is a contrast column:
 - '0' or '': the sample is not included in the analysis;
 - '1': the sample belongs to the **control** group;
 - '2': the sample belongs to the **treatment** test case group.
- For chipseq, N=NarrowPeaks and B=BroadPeaks
- Providing the pipeline with a design file (-d or --design):

```
rnaseq.py -d mydesign.tsv
```

Readset & Design File Validation



- Getting the Readset and Design files right is the most important task as a user.
- To avoid pipeline errors due to file issues, we have a file validator, `muggicValidator.py`, that can let you know if your Readset or Design files have issues.
- To validate your Readset file:


```
python $MUGQIC_PIPELINES_HOME/utils/muggicValidator.py -r myreadsetfile
```

- To validate your Design file:

```
python $MUGQIC_PIPELINES_HOME/utils/muggicValidator.py -d mydesignfile
```

GenPipes Overview

GenPipes Requirements:



GenPipes requires setting up 3 files:

- Configuration “ini” file: contains parameters used by tools in the pipeline
- Readset file: contains details about your samples
- Design file: contains details about sample comparisons. **NOT needed by all pipelines.**

Steps



- The Steps parameter tells GenPipes which steps to run.
- To know what the steps are, you can use `<pipeline>.py -h`
- To specify the steps: `rnaseq.py -h`

```
rnaseq.py -h
```


```
## full pipeline:  
rnaseq.py -s 1-25  
  
## specify steps 1 to 7 and 12  
rnaseq.py -s 1-7,12
```



- Steps:

- 1- picard_sam_to_fastq
 - 2- trimmomatic
 - 3- merge_trimmomatic_stats
 - 4- star
 - 5- picard_merge_sam_files
 - 6- picard_sort_sam
 - 7- picard_mark_duplicates
 - 8- picard_rna_metrics
 - 9- estimate_ribosomal_rna
 - 10- bam_hard_clip
 - 11- rnaseq
 - 12- wiggle
 - 13- raw_counts
 - 14- raw_counts_metrics
 - 15- cufflinks
 - 16- cuffmerge
 - 17- cuffquant
 - 18- cuffdiff
 - 19- cuffnorm
 - 20- fpkm_correlation_matrix
 - 21- gq_seq_utils_exploratory_analysis_rnaseq
 - 22- differential_expression
 - 23- differential_expression_goseq
 - 24- ihec_metrics
 - 25- verify_bam_id

Running GenPipes: Putting it all together



```
## create the command script  
rnaseq.py -c rnaseq.base.ini -r readset.txt -d designfile.txt -s 1-15 -j slurm > commands.txt
```



pipeline



Ini file



readset file



design file



steps



scheduler



script

```
## run the script  
bash commands.txt
```

```
## to check commands in queue:  
squeue -u $USER
```


Running GenPipes: Output Files, Logs and Errors

- As jobs finish, log files will be created in the “**job_output**” folder.
- If everything is Ok, you will get **MUGQICexitStatus:0** or **Exit_status=0**
- If you get an error that you cannot figure out, you can send us an email.
- For every analysis step, a folder will be created with the appropriate output files. Look into all the folders to know what each pipeline produces.
- If you made a mistake and would like to cancel ALL jobs in your queue, use:

```
scancel -u $USER
```

Running GenPipes:

Interpreting the output




- Our pipelines are a framework that runs third party software that is used in particular fields. We do not write the software ourselves. We control the framework that runs the tools.
- For knowledge about the format of the output, you need to read the documentation of the tools that produced the output.
- You can contact us for information that you can not find.

A note about your laptop and the server



- Your laptop has a graphic user interface and useful software like Excel but rarely has enough capacity to run a full bioinformatics analysis.
- The server has enough capacity to run analysis but has no graphic user interface and can be hard to create files, like the readset and design files, and to visualize plots.
- What usually happens is that we create and view files on our laptops but run larger computations on the server. So moving files between the server and your laptop is useful. In this workshop, we will use CyberDuck to transfer between the two systems.



A note about different Operating Systems

- If you transfer files (Readset or Design files) between Mac, Linux and Windows, sometimes you get errors related to differences in how files are handled by the different systems, specifically relating to end of file (EOF) characters.
- In order to fix file formats to suit the unix environment of most servers, try [dos2unix](#):

```
dos2unix myfile.txt
```

- Another thing to try in vim is:

```
:%s/\r/\r/g
```

Bioinformatics Resources



- Our website: <http://www.computationalgenomics.ca/>
- Our github pages: <http://c3g.github.io>
- Our email: support@computationalgenomics.ca
- C3G Open Door: <http://www.computationalgenomics.ca/open-door/>
- Guillimin Training: <http://www.hpc.mcgill.ca/index.php/training>
- Online Forums: BioStars, SeqAnswers, ...
- Various Workshops/Resources: C3G, bioinformatics.ca, Coursera courses, Software Carpentry ...

Questions?



Canadian Centre for
Computational
Genomics



Fear is your worst enemy... Don't be intimidated;
experiment, test, fail and learn!

Backup First though!



Genome Québec



Genome Canada



Ontario Genomics



McGill



Canadian Centre for
Computational
Genomics

SickKids®

MiC:M McGill initiative in
Computational Medicine

Thank you!

 HPC4Health

 **IHEC**
International Human Epigenome Consortium

compute | calcul
canada | canada


bioinformatics.ca

 Terry Fox
PROFYLE


GenAP





Calcul Québec


canarie

Data from Nanuq

- csv file from Nanuq:



Project **Reads** Sample Submission (1) Samples (0) Libraries (15) **HiSeq Read Sets (18)** Documents (0) Analyses (0)

CSV Download Read Files Help with icons BAM Files Help


Read Sets (18 elements)

Filter

<input type="checkbox"/>	Name	Alias	Barcode	Library Barcode	Quote(s)	Library Source	Library Type	Type of Sequencing	Run Concentration (pM)	Concentration Method
--------------------------	------	-------	---------	-----------------	----------	----------------	--------------	--------------------	------------------------	----------------------

- This file is very useful in helping you find your experimental data to prepare the readset file

Download your data from Nanuq



Project Details | Sample Submission (1) | Samples (0) | Libraries (15) | **HiSeq Read Sets (18)** | Documents (0) | Analyses (0)

CSV | **Download Read Files** | Help with icons | BAM Files Help

Read Sets (18 elements)

Filter

<input type="checkbox"/>	Name	Alias	Barcode	Library Barcode	Quote(s)	Library Source	Library Type	Type of Sequencing	Run Concentration (pM)	Concentration Method
--------------------------	------	-------	---------	-----------------	----------	----------------	--------------	--------------------	------------------------	----------------------

- if downloading all readsets in the project, copy and paste download command to the server



Download Type

Download all read files for the technology 

Download files from selected reads

Download all read files for the technology

To download all the read files of the current technology, copy and paste the Unix command below. Replace the "<username>" and "<password>" placeholders with your Nanuq credentials. The download will resume any interrupted download and will not restart completed files.


NOTE: Your username and password may be visible to other users of the system running the command.

```
wget -O - "https://genomequebec.mcgill.ca/nanuqMPS/readsetList?projectId=13340&tech=HiSeq" --no-cookies --no-check-certificate --post-data '_username=username&_password=password' | wget --no-cookies --no-check-certificate --post-data '_username=username&_password=password' &
```



Download your data from Nanuq

- if downloading selected readsets in the project:



Download Type

Download all read files for the technology

Download files from selected reads

Download files from selected reads

NOTE: If one or more files have to be retrieved from the tape archive, the download may take a few minutes longer to complete.

You may choose between 3 types of download:

- Compressed File:** contains all the files of the chosen type for the selected reads.
- Text File:** contains URL links to files of the chosen type for the selected reads including md5 files.
- Md5 File:** contains all the md5 signature of the chosen type for the selected reads.

Type of Download:

Compressed file

Text file with URL links (UNIX/Linux/Mac OS X users only)

Md5 File

Type of Files to Download: Fastq R1 Fastq R2 Bam Bam Index

Close Download

- Download the file, and follow the Readme.txt in order to download your files.